

Product datasheet for GA217090

Mouse Marchf4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	March4 CRISPRa kit - CRISPR gene activation of mouse membrane-associated ring finger (C3HC4) 4
Symbol:	Marchf4
Components:	GA217090G1, Marchf4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GEI00074) GA217090G2, Marchf4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GEI00074) GA217090G3, Marchf4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GEI00074) 1 CRISPRa-Enhancer vector, SKU GEI00056 1 CRISPRa scramble vector, SKU GEI00077
UniProt ID:	<u>Q80TE3</u>
Synonyms:	BC056494; mKIAA1399
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

