

Product datasheet for **GA217105**

Mouse Lonrf2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Lonrf2 CRISPRa kit - CRISPR gene activation of mouse LON peptidase N-terminal domain and ring finger 2
Symbol:	Lonrf2
Components:	GA217105G1, Lonrf2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217105G2, Lonrf2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217105G3, Lonrf2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	2900060P06Rik; AI851349
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)