


Product datasheet for **GA217113**

Mouse A530006G24Rik activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	A530006G24Rik CRISPRa kit - CRISPR gene activation of mouse RIKEN cDNA A530006G24 gene
Symbol:	A530006G24Rik
Components:	GA217113G1, A530006G24Rik gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217113G2, A530006G24Rik gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217113G3, A530006G24Rik gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)