

Product datasheet for **GA217137**

Mouse Aadacl4fm1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	9430007A20Rik CRISPRa kit - CRISPR gene activation of mouse RIKEN cDNA 9430007A20 gene
Symbol:	Aadacl4fm1
Components:	GA217137G1, Aadacl4fm1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217137G2, Aadacl4fm1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217137G3, Aadacl4fm1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BM81</u>
Synonyms:	Aadacl4; Gm1034
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)