

Product datasheet for **GA217144**

Mouse *Adgrf3* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Adgrf3</i> CRISPRa kit - CRISPR gene activation of mouse adhesion G protein-coupled receptor F3
Symbol:	<i>Adgrf3</i>
UniProt ID:	Q58Y75
Synonyms:	Gm1041; Gpr113; Pgr23
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)