


Product datasheet for **GA217146**

Mouse Gm1673 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Gm1673 CRISPRa kit - CRISPR gene activation of mouse predicted gene 1673 |
| Symbol: | Gm1673 |
| Components: | GA217146G1, Gm1673 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217146G2, Gm1673 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217146G3, Gm1673 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID: | <u>Q3UR78</u> |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online »](#)