


Product datasheet for **GA217152**

Mouse A330070K13Rik activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	A330070K13Rik CRISPRa kit - CRISPR gene activation of mouse RIKEN cDNA A330070K13 gene
Symbol:	A330070K13Rik
Components:	GA217152G1, A330070K13Rik gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217152G2, A330070K13Rik gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217152G3, A330070K13Rik gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	Gm1051
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)