

Product datasheet for **GA217161**

Mouse Pttg1ip2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	1700015F17Rik CRISPRa kit - CRISPR gene activation of mouse RIKEN cDNA 1700015F17 gene
Symbol:	Pttg1ip2
Components:	GA217161G1, Pttg1ip2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217161G2, Pttg1ip2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217161G3, Pttg1ip2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>D3YUK8</u>
Synonyms:	Gm1686
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)