

Product datasheet for GA217163

Mouse Lrrc43 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Lrrc43 CRISPRa kit - CRISPR gene activation of mouse leucine rich repeat containing 43
Symbol:	Lrrc43
UniProt ID:	Q3V0L5
Synonyms:	Gm1061
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)