

Product datasheet for **GA217168**

Mouse Igkv16-104 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Igkv16-104 CRISPRa kit - CRISPR gene activation of mouse immunoglobulin kappa variable 16-104
Symbol:	Igkv16-104
Components:	GA217168G1, Igkv16-104 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217168G2, Igkv16-104 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217168G3, Igkv16-104 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)