

Product datasheet for **GA217173**

Mouse 4930590J08Rik activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	4930590J08Rik CRISPRa kit - CRISPR gene activation of mouse RIKEN cDNA 4930590J08 gene
Symbol:	4930590J08Rik
Components:	GA217173G1, 4930590J08Rik gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217173G2, 4930590J08Rik gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217173G3, 4930590J08Rik gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8CDN1</u>
Synonyms:	Gm1071
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)