

Product datasheet for **GA217175**

Mouse Tsen2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tsen2 CRISPRa kit - CRISPR gene activation of mouse tRNA splicing endonuclease subunit 2
Symbol:	Tsen2
Components:	GA217175G1, Tsen2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217175G2, Tsen2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217175G3, Tsen2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6P7W5</u>
Synonyms:	AU067695
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)