

## Product datasheet for **GA217183**

### Mouse Lockd activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Lockd CRISPRa kit - CRISPR gene activation of mouse lncRNA downstream of Cdkn1b
Symbol:	Lockd
Components:	GA217183G1, Lockd gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217183G2, Lockd gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217183G3, Lockd gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)