

## Product datasheet for **GA217186**

### Mouse Sbk3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Sbk3 CRISPRa kit - CRISPR gene activation of mouse SH3 domain binding kinase family, member 3
Symbol:	Sbk3
Components:	GA217186G1, Sbk3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217186G2, Sbk3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217186G3, Sbk3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P0C5K0</u>
Synonyms:	Gm1078; Sgk110
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)