

Product datasheet for **GA217187**

Mouse Sbk2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Sbk2 CRISPRa kit - CRISPR gene activation of mouse SH3-binding domain kinase family, member 2
Symbol:	Sbk2
Components:	GA217187G1, Sbk2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217187G2, Sbk2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217187G3, Sbk2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P0C5K1</u>
Synonyms:	Gm1079; Sgk069
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)