

Product datasheet for **GA217189**

Mouse Rnf225 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Rnf225 CRISPRa kit - CRISPR gene activation of mouse ring finger protein 225
Symbol:	Rnf225
Components:	GA217189G1, Rnf225 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217189G2, Rnf225 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217189G3, Rnf225 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9D7D1</u>
Synonyms:	2310014L17Rik; AV237412; C77118
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)