

## Product datasheet for **GA217199**

### Mouse Lmtk3 activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | Lmtk3 CRISPRa kit - CRISPR gene activation of mouse lemur tyrosine kinase 3 |
| Symbol:       | Lmtk3   |
| UniProt ID:   | <a href="#">Q5XIV6</a>  |
| Synonyms:     | aatyk3; BC059845  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)     |

#### Product images:



[View online >](#)