


Product datasheet for **GA217217**

Mouse Ccdc33 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ccdc33 CRISPRa kit - CRISPR gene activation of mouse coiled-coil domain containing 33
Symbol:	Ccdc33
Components:	GA217217G1, Ccdc33 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217217G2, Ccdc33 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217217G3, Ccdc33 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q3ULW6
Synonyms:	4930535E21Rik; Gm1117; Tzip1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)