

Product datasheet for **GA217237**

Mouse Brwd3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Brwd3 CRISPRa kit - CRISPR gene activation of mouse bromodomain and WD repeat domain containing 3
Symbol:	Brwd3
Components:	GA217237G1, Brwd3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217237G2, Brwd3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217237G3, Brwd3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A2AHJ4
Synonyms:	Brod1; D030064D06Rik; Gm596
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)