

Product datasheet for **GA217238**

Mouse Cpxcr1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cpxcr1 CRISPRa kit - CRISPR gene activation of mouse CPX chromosome region, candidate 1
Symbol:	Cpxcr1
Components:	GA217238G1, Cpxcr1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217238G2, Cpxcr1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217238G3, Cpxcr1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3V0P1</u>
Synonyms:	Gm1143
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)