

Product datasheet for **GA217241**

Mouse Bclaf3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Bclaf3 CRISPRa kit - CRISPR gene activation of mouse Bclaf1 and Thrap3 family member 3
Symbol:	Bclaf3
Components:	GA217241G1, Bclaf3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217241G2, Bclaf3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217241G3, Bclaf3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A2AG58
Synonyms:	2310067E24Rik; Gm1146
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)