

Product datasheet for **GA217245**

Mouse Gm5169 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Gm5169 CRISPRa kit - CRISPR gene activation of mouse predicted gene 5169
Symbol:	Gm5169
Components:	GA217245G1, Gm5169 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217245G2, Gm5169 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217245G3, Gm5169 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	382277; Xmr
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)