


Product datasheet for **GA217250**

Mouse Gm5177 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Gm5177 CRISPRa kit - CRISPR gene activation of mouse glyceraldehyde-3-phosphate dehydrogenase pseudogene
Symbol:	Gm5177
Components:	GA217250G1, Gm5177 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217250G2, Gm5177 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217250G3, Gm5177 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)