

Product datasheet for **GA217251**

Mouse H2bu1-ps activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Hist3h2bb-ps CRISPRa kit - CRISPR gene activation of mouse histone cluster 3, H2bb, pseudogene
Symbol:	H2bu1-ps
Components:	GA217251G1, H2bu1-ps gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217251G2, H2bu1-ps gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217251G3, H2bu1-ps gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	4930534G10Rik; 4933432H21; B2Hist3h2bb; Hist3h2bb; R75370
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)