


Product datasheet for **GA217254**

Mouse Cd300ld3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cd300ld3 CRISPRa kit - CRISPR gene activation of mouse CD300 molecule like family member D3
Symbol:	Cd300ld3
Components:	GA217254G1, Cd300ld3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217254G2, Cd300ld3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217254G3, Cd300ld3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q6S1Q5
Synonyms:	AY457049; Cd300lh; clm-3; Clm3; Lmir7
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)