

Product datasheet for **GA217255**

Mouse Pfn4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pfn4 CRISPRa kit - CRISPR gene activation of mouse profilin family, member 4
Symbol:	Pfn4
UniProt ID:	Q9D6I3
Synonyms:	2900024P18Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)