


Product datasheet for **GA217268**

Mouse Tvp23a activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tvp23a CRISPRa kit - CRISPR gene activation of mouse trans-golgi network vesicle protein 23A
Symbol:	Tvp23a
Components:	GA217268G1, Tvp23a gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217268G2, Tvp23a gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217268G3, Tvp23a gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	Fam18a
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)