

Product datasheet for **GA217270**

Mouse Kctd16 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Kctd16 CRISPRa kit - CRISPR gene activation of mouse potassium channel tetramerisation domain containing 16
Symbol:	Kctd16
Components:	GA217270G1, Kctd16 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217270G2, Kctd16 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217270G3, Kctd16 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5DTY9</u>
Synonyms:	4930434H12Rik; Gm1267
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)