

## Product datasheet for **GA217279**

### Mouse PrameI5 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                  |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                     |
| Description:  | PrameI5 CRISPRa kit - CRISPR gene activation of mouse preferentially expressed antigen in melanoma like 5                                                                                                                                                                                                        |
| Symbol:       | PrameI5                                                                                                                                                                                                                                                                                                          |
| Components:   | GA217279G1, PrameI5 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA217279G2, PrameI5 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA217279G3, PrameI5 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| Synonyms:     | OTTMUSG00000010540                                                                                                                                                                                                                                                                                               |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                          |

#### Product images:



[View online »](#)