


Product datasheet for **GA217289**

Mouse Igkv8-30 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Igkv8-30 CRISPRa kit - CRISPR gene activation of mouse immunoglobulin kappa chain variable 8-30
Symbol:	Igkv8-30
Components:	GA217289G1, Igkv8-30 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217289G2, Igkv8-30 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217289G3, Igkv8-30 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)