


Product datasheet for **GA217293**

Mouse Ceacam3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ceacam3 CRISPRa kit - CRISPR gene activation of mouse carcinoembryonic antigen-related cell adhesion molecule 3
Symbol:	Ceacam3
Components:	GA217293G1, Ceacam3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217293G2, Ceacam3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217293G3, Ceacam3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	cea12; EG384557; Psg24
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)