

Product datasheet for **GA217296**

Mouse Scgb1b3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Scgb1b3 CRISPRa kit - CRISPR gene activation of mouse secretoglobin, family 1B, member 3
Symbol:	Scgb1b3
Components:	GA217296G1, Scgb1b3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217296G2, Scgb1b3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217296G3, Scgb1b3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	Abpa3; Abpa21; EG384585; Gm5325
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)