

Product datasheet for **GA217298**

Mouse Kash5 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ccdc155 CRISPRa kit - CRISPR gene activation of mouse coiled-coil domain containing 155
Symbol:	Kash5
Components:	GA217298G1, Kash5 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217298G2, Kash5 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217298G3, Kash5 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q80VJ8</u>
Synonyms:	Gm1434; KASH5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)