


Product datasheet for **GA217306**

Mouse Gm1943 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Gm1943 CRISPRa kit - CRISPR gene activation of mouse WD repeat domain 70 pseudogene
Symbol:	Gm1943
Components:	GA217306G1, Gm1943 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217306G2, Gm1943 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217306G3, Gm1943 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)