

Product datasheet for **GA217314**

Mouse Igkv4-54 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Igkv4-54 CRISPRa kit - CRISPR gene activation of mouse immunoglobulin kappa chain variable 4-54
Symbol:	Igkv4-54
Components:	GA217314G1, Igkv4-54 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217314G2, Igkv4-54 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217314G3, Igkv4-54 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)