

Product datasheet for **GA217315**

Mouse H2aI3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	H2aI3 CRISPRa kit - CRISPR gene activation of mouse H2A histone family member L3
Symbol:	H2aI3
Components:	GA217315G1, H2aI3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217315G2, H2aI3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217315G3, H2aI3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	4930557A04Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)