


Product datasheet for **GA217316**

Mouse H2a1n activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	H2a1n CRISPRa kit - CRISPR gene activation of mouse H2A histone family member L1N
Symbol:	H2a1n
Components:	GA217316G1, H2a1n gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217316G2, H2a1n gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217316G3, H2a1n gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	EG385328; Gm5382
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)