

Product datasheet for **GA217320**

Mouse **Tex28** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Tex28 CRISPRa kit - CRISPR gene activation of mouse testis expressed 28
Symbol:	Tex28
Components:	GA217320G1, Tex28 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217320G2, Tex28 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217320G3, Tex28 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	Gm1538; Gm1811
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)