

## Product datasheet for **GA217321**

### Mouse AU015836 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                     |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                        |
| Description:  | AU015836 CRISPRa kit - CRISPR gene activation of mouse expressed sequence AU015836                                                                                                                                                                                                                                  |
| Symbol:       | AU015836                                                                                                                                                                                                                                                                                                            |
| Components:   | GA217321G1, AU015836 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA217321G2, AU015836 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA217321G3, AU015836 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                             |

#### Product images:



[View online »](#)