

## Product datasheet for **GA217324**

### Mouse Lca5l activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                            |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                               |
| Description:  | Lca5l CRISPRa kit - CRISPR gene activation of mouse Leber congenital amaurosis 5-like                                                                                                                                                                                                                      |
| Symbol:       | Lca5l                                                                                                                                                                                                                                                                                                      |
| Components:   | GA217324G1, Lca5l gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA217324G2, Lca5l gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA217324G3, Lca5l gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q8C0X0</a></u>                                                                                                                                                                                                                                                                              |
| Synonyms:     | 4921526F01Rik; Lcca5l                                                                                                                                                                                                                                                                                      |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                    |

#### Product images:



[View online »](#)