

## Product datasheet for **GA217338**

### Mouse *Tas2r102* activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Tas2r102 CRISPRa kit - CRISPR gene activation of mouse taste receptor, type 2, member 102
Symbol:	Tas2r102
Components:	GA217338G1, Tas2r102 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217338G2, Tas2r102 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217338G3, Tas2r102 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	mGR02; mt2r51; Tas2r2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)