

Product datasheet for **GA217340**

Mouse **Tas2r110** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Tas2r110 CRISPRa kit - CRISPR gene activation of mouse taste receptor, type 2, member 110
Symbol:	Tas2r110
Components:	GA217340G1, Tas2r110 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217340G2, Tas2r110 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217340G3, Tas2r110 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q7M712
Synonyms:	mt2r57; T2R10; Tas2r10
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)