

## Product datasheet for **GA217341**

### Mouse **Tas2r118** activation kit by **CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	Tas2r118 CRISPRa kit - CRISPR gene activation of mouse taste receptor, type 2, member 118
Symbol:	Tas2r118
Components:	GA217341G1, Tas2r118 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217341G2, Tas2r118 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217341G3, Tas2r118 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P59529</u>
Synonyms:	Gm465; mGR18; mt2r40; T2R18; Tas2r16; Tas2r18
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)