

Product datasheet for **GA217366**

Mouse Tdpoz3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tdpoz3 CRISPRa kit - CRISPR gene activation of mouse TD and POZ domain containing 3
Symbol:	Tdpoz3
Components:	GA217366G1, Tdpoz3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217366G2, Tdpoz3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217366G3, Tdpoz3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q717B4</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)