

Product datasheet for **GA217371**

Mouse D030018L15Rik activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	D030018L15Rik CRISPRa kit - CRISPR gene activation of mouse nuclear receptor coactivator 2 pseudogene
Symbol:	D030018L15Rik
Components:	GA217371G1, D030018L15Rik gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217371G2, D030018L15Rik gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217371G3, D030018L15Rik gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)