

Product datasheet for GA217372

Mouse D030025E07Rik activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	D030025E07Rik CRISPRa kit - CRISPR gene activation of mouse RIKEN cDNA D030025E07 gene
Symbol:	D030025E07Rik
Components:	GA217372G1, D030025E07Rik gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217372G2, D030025E07Rik gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217372G3, D030025E07Rik gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)