


Product datasheet for **GA217376**

Mouse Tigd4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tigd4 CRISPRa kit - CRISPR gene activation of mouse tigger transposable element derived 4
Symbol:	Tigd4
Components:	GA217376G1, Tigd4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217376G2, Tigd4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217376G3, Tigd4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BUZ3</u>
Synonyms:	C130063O11Rik; Gm418
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)