


Product datasheet for **GA217377**

Mouse Plcxd1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Plcxd1 CRISPRa kit - CRISPR gene activation of mouse phosphatidylinositol-specific phospholipase C, X domain containing 1
Symbol:	Plcxd1
Components:	GA217377G1, Plcxd1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217377G2, Plcxd1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217377G3, Plcxd1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	A330045H12Rik; Gm576
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)