


Product datasheet for **GA217413**

Mouse *Iqgap3* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Iqgap3</i> CRISPRa kit - CRISPR gene activation of mouse IQ motif containing GTPase activating protein 3
Symbol:	<i>Iqgap3</i>
Synonyms:	AI593484; D030034H08
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)