

Product datasheet for **GA217440**

Mouse Znr3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Znr3 CRISPRa kit - CRISPR gene activation of mouse zinc and ring finger 3
Symbol:	Znr3
UniProt ID:	<u>Q5SSZ7</u>
Synonyms:	Gm1167
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)